

Product Datasheet

KLKB1 (H chain, Cleaved-Arg390) antibody (orb126592)

Description Rabbit Polyclonal to KLKB1 (H chain, Cleaved-Arg390).

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen Internal

Concentration 1mg/ml


Note For research use only

Clonality Polyclonal

Uniprot ID **P03952**

Dilution Range WB: 1:500-1:3000, ELISA: 1:10000

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts from H...