

Product Datasheet

Dipeptidyl-peptidase 1 (H chain, Cleaved-Arg394) antibody (orb126593)

Description Rabbit Polyclonal to Dipeptidyl-peptidase 1 (H chain, Cleaved-Arg394).

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications WB

Immunogen C-terminal

Concentration 1mg/ml


Note For research use only

Clonality Polyclonal

Uniprot ID **P53634**

Dilution Range WB: 1:500-1:3000, ELISA: 1:20000

Expiration Date 12 months from date of receipt.



Western blot analysis of extracts from J...