


Product Datasheet

KATNA1 Antibody (orb1265932)

Description	KATNA1 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	Western blot analysis in 293 cell line l...
Immunogen	This KATNA1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 59-85 amino acids from the N-terminal region of human KATNA1.	
Target	KATNA1	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	56 kDa	
Uniprot ID	O75449	
NCBI	O75449	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	