

Product Datasheet

Mycb Antibody (orb1265950)

Description

Mycb Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This Rat Mycb antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 120-146 amino acids from the C-terminal region of rat Mycb.

Target

Mycb

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

18 kDa

Uniprot ID

P15063

NCBI


P15063

Dilution Range

For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis in mouse bladder t...