


Product Datasheet

Avpr1b Antibody (orb1265951)

Description	Avpr1b Antibody	
Species/Host	Rabbit	
Reactivity	Mouse	Western blot analysis in mouse liver tis...
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This Mouse Avpr1b antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 216-242 amino acids from the Central region of mouse Avpr1b.	
Target	Avpr1b	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	47 kDa	
Uniprot ID	Q9WU02	
NCBI	Q9WU02	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	