


Product Datasheet

Cst3 Antibody (orb1265952)

Antibody Description Description		Cst3 Antibody 
Species/Host	Rabbit	Western blot analysis in mouse spleen ti...
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This Mouse Cst3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 39-65 amino acids from the Central region of mouse Cst3.	
Target	Cst3	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	16 kDa	
Uniprot ID	P21460	
NCBI	P21460	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	