


Product Datasheet

Caspase 3 (Cleaved-Ser29) antibody (orb126597)

Description	Rabbit polyclonal antibody to Caspase 3 (Cleaved-Ser29)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	N-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P42574
Dilution Range	WB: 1:500-1:3000, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Western
blot
analysis of
extracts
from C...