



Product Datasheet

FBXO8 Antibody (orb1265971)

Description	FBXO8 Antibody	
Species/Host	Rabbit	
Reactivity	Mouse	
Conjugation	Unconjugated	
Tested Applications	WB	Western blot analysis in mouse bladder t...
Immunogen	This FBXO8 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 168-195 amino acids from the Central region of human FBXO8.	
Target	FBXO8	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	37 kDa	
Uniprot ID	Q9NRD0	
NCBI	Q9NRD0	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	