

Product Datasheet

Caspase 3 (p17, Cleaved-Asp175) antibody (orb126609)

Description Rabbit polyclonal antibody to Caspase 3 (p17, Cleaved-Asp175)

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Internal

Concentration 1mg/ml


Note For research use only

Clonality Polyclonal

Uniprot ID **P42574**

Dilution Range WB: 1:500-1:3000, ELISA: 1:20000

Expiration Date 12 months from date of receipt.



Western
blot
analysis of
extracts
from C...