

Product Datasheet

PGLYRP4 Antibody (orb1266121)

Description

PGLYRP4 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This PGLYRP4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 242-268 amino acids from the Central region of human PGLYRP4.

Target

PGLYRP4

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

41 kDa

Uniprot ID

[Q96LB8](#)

NCBI


[Q96LB8](#)

Dilution Range


For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western Blot
at 1:1000
dilution +
HACAT ...



Western blot
analysis in
A549 cell line
...