


Product Datasheet

ABHD1 Antibody (orb1266130)

Description	ABHD1 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This ABHD1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 289-315 amino acids from the C-terminal region of human ABHD1.	Western blot analysis in MDA-MB468 cell ...
Target	ABHD1	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	45 kDa	
Uniprot ID	Q96SE0	
NCBI	Q96SE0	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	