



Product Datasheet

PCDHGA5 Antibody (orb1266214)

Description	PCDHGA5 Antibody	MDA-MB453 
Species/Host	Rabbit	130 95 72  Western blot analysis in MDA-MB453 cell ...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This PCDHGA5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 254-280 amino acids from the N-terminal region of human PCDHGA5.	
Target	PCDHGA5	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	101 kDa	
Uniprot ID	Q9Y5G8	
NCBI	Q9Y5G8	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	