

## Product Datasheet

### METTL6 Antibody (orb1266225)

## Description

METTL6 Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse

### Conjugation

Unconjugated

### Tested

WB

### Applications

### Immunogen

This METTL6 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 168-195 amino acids from the Central region of human METTL6.

### Target

METTL6

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

33 kDa

### Uniprot ID

[Q8TCB7](#)

### NCBI

[Q8TCB7](#)

### Dilution Range

For WB starting dilution is: 1:1000

### Expiration Date

12 months from date of receipt.

m.brain  
95

