


Product Datasheet

HAUS8 Antibody (orb1266290)

Description	HAUS8 Antibody	 Western blot analysis in A549 cell line ...
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This HAUS8 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 28-54 amino acids from the N-terminal region of human HAUS8.	
Target	HAUS8	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	45 kDa	
Uniprot ID	Q9BT25	
NCBI	Q9BT25	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	