


Product Datasheet

Gamma-glutamyltransferase 4 (H chain, Cleaved-Thr472) antibody (orb126632)

Description	Rabbit Polyclonal to Gamma-glutamyltransferase 4 (H chain, Cleaved-Thr472).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9UJ14
Dilution Range	WB: 1:500-1:3000, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from J...