


Product Datasheet

ENTPD4 Antibody (orb1266356)

Western Blot analysis. Description		ENTPD4 Antibody	
Species/Host	Rabbit		
Reactivity	Human, Mouse		
Conjugation	Unconjugated		
Tested Applications	WB		Western blot analysis in mouse spleen ti...
Immunogen	This ENTPD4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 40-69 amino acids from the N-terminal region of human ENTPD4.		
Target	ENTPD4		
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.		
Form/Appearance	Liquid		
Concentration	batch dependent		
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.		
Note	For research use only		
Application notes	For WB starting dilution is: 1:1000		
Isotype	Rabbit Ig		
Clonality	Polyclonal		
MW	70 kDa		
Uniprot ID	Q9Y227		
NCBI	Q9Y227		
Dilution Range	For WB starting dilution is: 1:1000		
Expiration Date	12 months from date of receipt.		