

Product Datasheet

TFAP2D Antibody (orb1266524)

Description

TFAP2D Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, WB

Applications

Immunogen

This TFAP2D antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 308-336 amino acids from the C-terminal region of human TFAP2D.

Target

TFAP2D

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IF starting dilution is: 1:10~50For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

50 kDa

Uniprot ID

[Q7Z6R9](#)

NCBI



[Q7Z6R9](#)

Dilution Range


For IF starting dilution is: 1:10~50For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Fluorescent image of A549 cell stained w...



Western blot analysis in HL-60 cell line...