



Product Datasheet

THEM5 Antibody (orb1266577)

Description	THEM5 Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	Western blot analysis in A375 cell line ...
Immunogen	This THEM5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 71-99 amino acids from the Central region of human THEM5.	
Target	THEM5	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	28 kDa	
Uniprot ID	Q8N1Q8	
NCBI	Q8N1Q8	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	