


## Product Datasheet

### RNASE11 Antibody (orb1266580)

|                                                          |                                                                                                                                                                                                                        |                  |                                                                                                                                                                                                                  |
|----------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <div>Western Blot analysis.</div> <div>Description</div> |                                                                                                                                                                                                                        | RNASE11 Antibody | <div><div>MCF-7</div></div>                                                                                                   |
| Species/Host                                             | Rabbit                                                                                                                                                                                                                 |                  | <div><div>20</div><div>17</div><div>10</div><div>157982101</div></div> <div>Western blot analysis in MCF-7 cell line...</div> |
| Reactivity                                               | Human                                                                                                                                                                                                                  |                  |                                                                                                                                                                                                                  |
| Conjugation                                              | Unconjugated                                                                                                                                                                                                           |                  |                                                                                                                                                                                                                  |
| Tested Applications                                      | WB                                                                                                                                                                                                                     |                  |                                                                                                                                                                                                                  |
| Immunogen                                                | This RNASE11 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 68-95 amino acids from the Central region of human RNASE11.                                                  |                  |                                                                                                                                                                                                                  |
| Target                                                   | RNASE11                                                                                                                                                                                                                |                  |                                                                                                                                                                                                                  |
| Preservatives                                            | Supplied in PBS with 0.09% (W/V) sodium azide.                                                                                                                                                                         |                  |                                                                                                                                                                                                                  |
| Form/Appearance                                          | Liquid                                                                                                                                                                                                                 |                  |                                                                                                                                                                                                                  |
| Concentration                                            | batch dependent                                                                                                                                                                                                        |                  |                                                                                                                                                                                                                  |
| Storage                                                  | Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures. |                  |                                                                                                                                                                                                                  |
| Note                                                     | For research use only                                                                                                                                                                                                  |                  |                                                                                                                                                                                                                  |
| Application notes                                        | For WB starting dilution is: 1:1000                                                                                                                                                                                    |                  |                                                                                                                                                                                                                  |
| Isotype                                                  | Rabbit Ig                                                                                                                                                                                                              |                  |                                                                                                                                                                                                                  |
| Clonality                                                | Polyclonal                                                                                                                                                                                                             |                  |                                                                                                                                                                                                                  |
| MW                                                       | 22 kDa                                                                                                                                                                                                                 |                  |                                                                                                                                                                                                                  |
| Uniprot ID                                               | Q8TAA1                                                                                                                                                                                                                 |                  |                                                                                                                                                                                                                  |
| NCBI                                                     | Q8TAA1                                                                                                                                                                                                                 |                  |                                                                                                                                                                                                                  |
| Dilution Range                                           | For WB starting dilution is: 1:1000                                                                                                                                                                                    |                  |                                                                                                                                                                                                                  |
| Expiration Date                                          | 12 months from date of receipt.                                                                                                                                                                                        |                  |                                                                                                                                                                                                                  |