

Product Datasheet

DNAJC5 Antibody (orb1266692)

Description

DNAJC5 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This DNAJC5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 52-80 amino acids from the Central region of human DNAJC5.

Target

DNAJC5

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

22 kDa

Uniprot ID

[Q9H3Z4](#)

NCBI


[Q9H3Z4](#)

Dilution Range


For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis in Y79 cell line l...



Western blot analysis in mouse NIH-3T3 c...