

Product Datasheet

SNTB1 Antibody (orb1266698)

Description

SNTB1 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This SNTB1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 243-271 amino acids from the Central region of human SNTB1.

Target

SNTB1

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

58 kDa

Uniprot ID

[Q13884](#)

NCBI

[Q13884](#)


Dilution Range

For WB starting dilution is: 1:1000


Expiration Date

12 months from date of receipt.

NCI-H460



Western blot analysis in NCI-H460 cell l...



Western blot analysis in mouse spleen ti...