

Product Datasheet

TAS2R9 Antibody (orb1266922)

Description

TAS2R9 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This TAS2R9 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 154-180 amino acids from the Central region of human TAS2R9.

Target

TAS2R9

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

36 kDa

Uniprot ID

[Q9NYW1](#)

NCBI

[Q9NYW1](#)


Dilution Range

For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.

MDA-MB453



Western blot analysis in MDA-MB453 cell ...