

Product Datasheet

KLHL12 Antibody (orb1266934)

Description

KLHL12 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This KLHL12 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 376-404 amino acids from the C-terminal region of human KLHL12.

Target

KLHL12

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

63 kDa


Uniprot ID
[Q53G59](#)
NCBI
[Q53G59](#)
Dilution Range

For WB starting dilution is: 1:1000


Expiration Date

12 months from date of receipt.

m.NIH-3T3



Western blot analysis in mouse NIH-3T3 c...



Western blot analysis in HeLa cell line ...