

Product Datasheet

Stk11 Antibody (orb1266943)

Description

Stk11 Antibody

Species/Host

Rabbit

Reactivity

Mouse

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

This Mouse Stk11 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 393-420 amino acids from the C-terminal region of mouse Stk11.

Target

Stk11

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For IF starting dilution is: 1:25For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:500-1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

49 kDa

Uniprot ID

[Q9WTK7](#)

NCBI


[Q9WTK7](#)

Dilution Range


For IF starting dilution is: 1:25For IHC-P starting dilution is: 1:25For WB starting dilution is: 1:500-1:1000

Expiration Date


12 months from date of receipt.



Fluorescent image of A549 cells stained ...



Immunohistochemical analysis of paraffin...



Western blot analysis of lysate from NCI...