

## Product Datasheet

### PLEKHG5 Antibody (orb1266963)

**Description**

PLEKHG5 Antibody

**Species/Host**

Rabbit

**Reactivity**

Human

**Conjugation**

Unconjugated

**Tested**

WB

**Applications**
**Immunogen**

This PLEKHG5 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 175-202 amino acids from the N-terminal region of human PLEKHG5.

**Target**

PLEKHG5

**Preservatives**

Supplied in PBS with 0.09% (W/V) sodium azide.

**Form/Appearance**

Liquid

**Concentration**

batch dependent

**Storage**

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

**Note**

For research use only

**Application notes**

For WB starting dilution is: 1:1000

**Isotype**

Rabbit Ig

**Clonality**

Polyclonal

**MW**

117 kDa

**Uniprot ID**

[O94827](#)

**NCBI**


[O94827](#)

**Dilution Range**

For WB starting dilution is: 1:1000

**Expiration Date**

12 months from date of receipt.



Western blot analysis in A549 cell line ...