

Product Datasheet

TSPY3 Antibody (orb1266964)

Description

TSPY3 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This TSPY3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 277-305 amino acids from the C-terminal region of human TSPY3.

Target

TSPY3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


35 kDa

Uniprot ID
P0CV98
NCBI
P0CV98
Dilution Range


For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western Blot
at 1:1000
dilution +
human ...



Western blot
analysis in
293 cell line
l...