




Product Datasheet

IL17D Antibody (orb1267121)

Description	IL17D Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	FC, IHC-P, WB	Western blot analysis of IL27 using rabb...
Immunogen	This IL-27 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 57-87 amino acids from the Central region of human IL-27.	
Target	IL17D	IL-27 Antibody immunohistochemistry anal...
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	Flow cytometric analysis of 293 cells (r...
Clonality	Polyclonal	
MW	22 kDa	
Uniprot ID	Q8TAD2	
NCBI	Q8TAD2	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	