



Product Datasheet

PTP4A1 Antibody (orb1267241)

Description	PTP4A1 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	WB	Western blot analysis in mouse bladder t...
Immunogen	This PTP4A1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 140-168 amino acids from the C-terminal region of human PTP4A1.	
Target	PTP4A1	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	Western blot analysis of PTP4A1 using ra...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	20 kDa	
Uniprot ID	Q93096	
NCBI	Q93096	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	