



Product Datasheet

CLCN2 Antibody (orb1267370)

Description	CLCN2 Antibody	
Species/Host	Rabbit	
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This CLCN2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 209-237 amino acids from the N-terminal region of human CLCN2.	 Western blot analysis in HL-60 cell line...
Target	CLCN2	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	99 kDa	
Uniprot ID	P51788	
NCBI	P51788	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	 Western blot analysis in mouse Neuro-2a ...