

Product Datasheet

H2BFWT Antibody (orb1267416)

Description

H2BFWT Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This H2BFWT antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 62-88 amino acids from the Central region of human H2BFWT.

Target

H2BFWT

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

20 kDa

Uniprot ID

[Q7Z2G1](#)

NCBI


[Q7Z2G1](#)

Dilution Range

For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western Blot
at 1:1000
dilution +
human ...