

Product Datasheet

TAS2R31 Antibody (orb1267455)

Description

TAS2R31 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications

Immunogen

This TAS2R31 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 280-309 amino acids from the C-terminal region of human TAS2R31.

Target

TAS2R31

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

35 kDa

Uniprot ID

P59538

NCBI


P59538

Dilution Range

For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western Blot
at 1:1000
dilution +
MDA-MB...