

Product Datasheet

TSC22D2 Antibody (orb1267500)

Description

TSC22D2 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, WB

Applications

Immunogen

This TSC22D2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 34-62 amino acids from the N-terminal region of human TSC22D2.

Target

TSC22D2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IF starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

79 kDa

Uniprot ID

[O75157](#)

NCBI


[O75157](#)

Dilution Range


For WB starting dilution is: 1:1000For IF starting dilution is: 1:10~50

Expiration Date

12 months from date of receipt.



Western blot analysis in HepG2 cell line...



Fluorescent confocal image of A431 cell ...