

## Product Datasheet

### TMPRSS11E Antibody (orb1267505)

## Description

TMPRSS11E Antibody

### Species/Host

Rabbit

### Reactivity

Human

### Conjugation

Unconjugated

### Tested

WB

### Applications

### Immunogen

This TMPRSS11E antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 249-277 amino acids from the Central region of human TMPRSS11E.

### Target

TMPRSS11E

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

48 kDa

### Uniprot ID

[Q9UL52](#)

### NCBI


[Q9UL52](#)

### Dilution Range

For WB starting dilution is: 1:1000

### Expiration Date

12 months from date of receipt.



Western blot analysis in HepG2 cell line...