



Product Datasheet

POLD3 Antibody (orb1267612)

Description	POLD3 Antibody	
Species/Host	Rabbit	 Western blot analysis in Jurkat cell lin...
Reactivity	Human, Mouse	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This POLD3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 416-443 amino acids from the C-terminal region of human POLD3.	 Western blot analysis in mouse spleen ti...
Target	POLD3	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	51 kDa	
Uniprot ID	Q15054	
NCBI	Q15054	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	