

Product Datasheet

SIGLEC12 Antibody (orb1267623)

Description

SIGLEC12 Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications

Immunogen

This SIGLEC12 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 565-595 amino acids from the C-terminal region of human SIGLEC12.

Target

SIGLEC12

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

65 kDa

Uniprot ID

[Q96PQ1](#)

NCBI

[Q96PQ1](#)


Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100


Expiration Date

12 months from date of receipt.

m.bladder
130



Western blot analysis of SIGLEC12 Antibody...



Formalin-fixed and paraffin-embedded hum...