


Product Datasheet

DLEC1 Antibody (orb1267666)

Antibody Bioscients. Description		DLEC1 Antibody	 Western blot analysis in NCI-H460 cell l...
Species/Host		Rabbit	
Reactivity		Human	
Conjugation		Unconjugated	
Tested Applications		WB	
Immunogen		This DLEC1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 263-292 amino acids from the N-terminal region of human DLEC1.	
Target		DLEC1	
Preservatives		Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance		Liquid	
Concentration		batch dependent	
Storage		Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note		For research use only	
Application notes		For WB starting dilution is: 1:1000	
Isotype		Rabbit Ig	
Clonality		Polyclonal	
MW		196 kDa	
Uniprot ID		Q9Y238	
NCBI		Q9Y238	
Dilution Range		For WB starting dilution is: 1:1000	
Expiration Date		12 months from date of receipt.	