


Product Datasheet

LINGO3 Antibody (orb1267737)

Description	LINGO3 Antibody	NCI-H460
Species/Host	Rabbit	 Western blot analysis in NCI-H460 cell l...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	WB	
Immunogen	This LINGO3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 395-424 amino acids from the C-terminal region of human LINGO3.	
Target	LINGO3	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	65 kDa	
Uniprot ID	P0C6S8	
NCBI	P0C6S8	
Dilution Range	For WB starting dilution is: 1:1000	
Expiration Date	12 months from date of receipt.	