


Product Datasheet

CAD (phospho-Thr456) antibody (orb126774)

Description	Rabbit Polyclonal to CAD (phospho-Thr456).
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Internal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P27708
Dilution Range	IHC-P: 1:50-1:100, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...