

Product Datasheet

Stk3 Antibody (orb1267788)

Description

Stk3 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This Mouse Stk3 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 387-415 amino acids from the C-terminal region of mouse Stk3.

Target

Stk3

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


57 kDa

Uniprot ID
[Q9JI10](https://www.uniprot.org/entry/Q9JI10)
NCBI
[Q9JI10](https://www.ncbi.nlm.nih.gov/nuclot/130319102)
Dilution Range


For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis in HeLa cell line ...



Western blot analysis of Stk3 using rabb...