



Product Datasheet

CBLN4 antibody (orb126779)

Description	Rabbit Polyclonal to CBLN4.
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q9NTU7
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts from C...



Immunohistochemical analysis of formalin...