


Product Datasheet

CTBP1 Antibody (orb1267852)


Description

CTBP1 Antibody


Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	This CTBP1 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 413-440 amino acids from the C-terminal region of human CTBP1.
Target	CTBP1
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.
Form/Appearance	Liquid
Concentration	batch dependent
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.
Note	For research use only
Application notes	For WB starting dilution is: 1:1000For IF starting dilution is: 1:10~50
Isotype	Rabbit Ig
Clonality	Polyclonal
MW	48 kDa
Uniprot ID	Q13363
NCBI	Q13363
Dilution Range	For WB starting dilution is: 1:1000For IF starting dilution is: 1:10~50
Expiration Date	12 months from date of receipt.



Western Blot
at 1:2000
dilution Lane
1: ...



Western blot
analysis in
Hela, HL-60,
Ju...



Fluorescent
confocal
image of
Hela cell ...