

Product Datasheet

TLR7 Antibody (orb1267880)

Description

TLR7 Antibody

150

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

IHC-P, WB

Applications
Immunogen

This TLR7 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 98-128 amino acids from the N-terminal region of human TLR7.

Target

TLR7

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

121 kDa

Uniprot ID
Q9NYK1
NCBI
Q9NYK1
Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:50~100

Expiration Date


12 months from date of receipt.



Antibody is used in Western blot to detect...



Western blot analysis in Ramos cell line...



Formalin-fixed and paraffin-embedded human tissue...