

Product Datasheet

POLG Antibody (orb1267926)

Description

POLG Antibody

Species/Host

Rabbit

Reactivity

Human, Mouse

Conjugation

Unconjugated

Tested

FC, IF, IHC-P, WB

Applications
Immunogen

This POLG antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1120-1148 amino acids from the C-terminal region of human POLG.

Target

POLG

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IF starting dilution is: 1:10~50For IHC-P starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

140 kDa

Uniprot ID

P54098

NCBI

P54098


Dilution Range

For WB starting dilution is: 1:1000For FACS starting dilution is: 1:10~50For IF starting dilution is: 1:10~50For IHC-P starting dilution is: 1:10~50


Expiration Date

12 months from date of receipt.


m.cerebellum
250



Western blot analysis in mouse cerebellu...



Flow cytometric analysis of CEM cells (r...



Confocal immunofluorescent analysis of P...