

Product Datasheet

CDH10 Antibody (orb1268021)

Description

CDH10 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

FC, IHC-P, WB

Applications

Immunogen

This CDH10 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 25-55 amino acids from the N-terminal region of human CDH10.

Target

CDH10

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

88 kDa

Uniprot ID

[Q9Y6N8](#)

NCBI


[Q9Y6N8](#)

Dilution Range

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50

Expiration Date


12 months from date of receipt.



Western blot analysis of CDH10 Antibody ...



Formalin-fixed and paraffin-embedded hum...



Flow cytometric analysis of CEM cells us...