

Product Datasheet

ALPI Antibody (orb1268035)

Description

ALPI Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

IF, IHC-P, WB

Applications

Immunogen

This Alkaline Phosphatase (ALPI) antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 274-304 amino acids from the Central region of human Alkaline Phosphatase (ALPI).

Target

ALPI

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:100For IF starting dilution is: 1:100

Isotype

Rabbit Ig

Clonality

Polyclonal

MW

57 kDa

Uniprot ID

[P09923](#)

NCBI


[P09923](#)

Dilution Range


For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:100For IF starting dilution is: 1:100

Expiration Date


12 months from date of receipt.



Western blot analysis of lysate from hum...



Immunohistochemical analysis of paraffin...



Western blot analysis of lysate from hum...