



Product Datasheet

MAP3K4 Antibody (orb1268053)

Description	MAP3K4 Antibody	
Species/Host	Rabbit	Western blot analysis in human normal Ut...
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	IHC-P, WB	
Immunogen	This MEKK4 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 1067-1097 amino acids from the Central region of human MEKK4.	
Target	MAP3K4	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	MEKK4 Antibody (Center C1081) immunohist...
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	182 kDa	
Uniprot ID	Q9Y6R4	
NCBI	Q9Y6R4	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50	
Expiration Date	12 months from date of receipt.	