



Product Datasheet

LXN Antibody (orb1268131)


Description	LXN Antibody	
Species/Host	Rabbit	
Reactivity	Human	
Conjugation	Unconjugated	
Tested Applications	FC, IHC-P, WB	
Immunogen	This Latexin antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 78-107 amino acids from the Central region of human Latexin.	Western Blot at 1:1000 dilution + MCF-7 ...
Target	LXN	
Preservatives	Supplied in PBS with 0.09% (W/V) sodium azide.	
Form/Appearance	Liquid	
Concentration	batch dependent	
Storage	Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.	Western blot analysis of Latexin using r...
Note	For research use only	
Application notes	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	
Isotype	Rabbit Ig	
Clonality	Polyclonal	
MW	26 kDa	
Uniprot ID	Q9BS40	
NCBI	Q9BS40	
Dilution Range	For WB starting dilution is: 1:1000For IHC-P starting dilution is: 1:10~50For FACS starting dilution is: 1:10~50	Formalin-fixed and paraffin-embedded hum...
Expiration Date	12 months from date of receipt.	



Western Blot at 1:1000 dilution + MCF-7 ...



Western blot analysis of Latexin using r...



Formalin-fixed and paraffin-embedded hum...