

## Product Datasheet

### IDE Antibody (orb1268222)

## Description

IDE Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse

### Conjugation

Unconjugated

### Tested

IHC-P, WB

### Applications

### Immunogen

This IDE antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 406-435 amino acids from the Central region of human IDE.

### Target

IDE

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:2000 For IHC-P starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

118 kDa

### Uniprot ID

[P14735](#)

### NCBI


[P14735](#)

### Dilution Range


For WB starting dilution is: 1:2000 For IHC-P starting dilution is: 1:10~50

### Expiration Date


12 months from date of receipt.



Western Blot  
at 1:2000  
dilution Lane  
1: ...



Western Blot  
at 1:2000  
dilution Lane  
1: ...



Western blot  
analysis of  
anti-IDE  
Antibo...