

Product Datasheet

NMT2 Antibody (orb1268224)

Description
NMT2 Antibody
Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

WB

Applications
Immunogen

This NMT2 antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 125-153 amino acids from the C-terminal region of human NMT2.

Target

NMT2

Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

Form/Appearance

Liquid

Concentration

batch dependent

Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

Note

For research use only

Application notes

For WB starting dilution is: 1:1000

Isotype

Rabbit Ig

Clonality

Polyclonal

MW


57 kDa

Uniprot ID
O60551
NCBI
O60551
Dilution Range


For WB starting dilution is: 1:1000

Expiration Date

12 months from date of receipt.



Western blot analysis in K562 cell line ...



Western blot analysis of NMT2 using rabb...