

## Product Datasheet

### IL9R Antibody (orb1268363)

## Description

IL9R Antibody

### Species/Host

Rabbit

### Reactivity

Human, Mouse

### Conjugation

Unconjugated

### Tested

IHC-P, WB

### Applications

### Immunogen

This IL9R antibody is generated from rabbits immunized with a KLH conjugated synthetic peptide between 170-199 amino acids from the Central region of human IL9R.

### Target

IL9R

### Preservatives

Supplied in PBS with 0.09% (W/V) sodium azide.

### Form/Appearance

Liquid

### Concentration

batch dependent

### Storage

Store at 4°C for three months and -20°C, stable for up to one year. As with all antibodies care should be taken to avoid repeated freeze thaw cycles. Antibodies should not be exposed to prolonged high temperatures.

### Note

For research use only

### Application notes

For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

### Isotype

Rabbit Ig

### Clonality

Polyclonal

### MW

57 kDa

### Uniprot ID

[Q01113](#)

### NCBI


[Q01113](#)

### Dilution Range


For WB starting dilution is: 1:1000 For IHC-P starting dilution is: 1:10~50

### Expiration Date

12 months from date of receipt.



Western Blot at  
1:1000 dilution +  
Jurkat...



IL9R Antibody  
immunohistochemistry  
analy...